VARISAN[®] A.T.E.[®]

VARISAN[®] RECOVER24

VARISAN[®] RECOVER35

- IT Istruzioni per l'uso
- EN Instructions for use
- FR Instructions d'utilisation
- ES Instrucciones de uso
- DE Gebrauchsanweisung
- PT Instruções de utilização
- PL Instrukcje użytkowania
- DA Instruktionser til brug
- NO Instruksjoner for bruk
- RU Инструкции по применению
- NL Gebruiksaanwijzing
- HR Upute za upotrebu
- EL Οδηγίες χρήσης
- AR تعليمات الاستخدام
- HI उपयोग के लिए निर्देश
- TR Kullanım talimatları

SCIENCE & MEDICAL
QUALITY

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EN

GENERAL INDICATIONS TO PUT ON THE STOCKINGS CORRECTLY

- Before putting on compression stockings, make sure that legs are clean and dry to prevent skin irritations and to make them easy to wear
- It is recommended to use rubber gloves to prevent rings or unfilled nails from damaging the stocking
- Turn the stocking inside out up to the heel
- Insert the foot, place the heel in the appropriate area of the stocking and make it slide along the leg without overstretching, accompanying the fabric with small upward movements

- Be sure that, once the stocking is on, the fabric does not form uncomfortable creases

INTENDED USE

VARISAN® A.T.E.* 18mmHg

Circular knit compression stockings for hospital use, suitable for the compression of lower limbs, mainly for the anti-thromboembolic preventive treatment.

VARISAN® RECOVER24 24mmHg

VARISAN® RECOVER35 35mmHg

Circular knit compression stockings for hospital use, suitable for the compression of lower limbs, mainly for the post-operation treatment of vascular pathologies or for the treatment of post-operation oedemas and haematomas.

TARGET PATIENT

The target patient may use the compression stockings after consulting with the attending doctor or professional in the sector and by following the supplier's instructions (measures, sizing, compression class).

INDICATIONS

VARISAN® A.T.E.* 18mmHg

- Treatment of pre-, intra-, and post-operative thrombosis of immobilized/bedridden patients
- Prevention and treatment of post-operative and post-traumatic oedema

VARISAN® RECOVER24 24mmHg

VARISAN® RECOVER35 35mmHg

- Short-term post-operative compression in vascular surgery, to promote venous return, and in post-operative compression in orthopaedic surgery
 - Compression treatment after acute superficial and deep vein thrombosis
 - Prevention and treatment of venous ulcers in mobilised/partially mobilised patients
 - Treatment of mild to moderate phlebo-lymphoedemas
- ## CONTRAINDICATIONS
- Peripheral artery occlusive disease with ankle-brachial index (ABI) below 0,7 or peripheral blood pressure below 80 mmHg or uncompensated
 - Acute skin infections
 - Exudative dermatitis
 - Acute dermo-hypodermatitis
 - Diseases of lower limbs
 - Oedemas of non-venous (dysmetabolic, diabetic, renal, cardiac, endocrine...) origin, in the acute phase or undiagnosed
 - Known history of contact allergy to medical compression stockings and their components
 - Congenital/acquired anatomical abnormalities (except for the possible use of custom-made braces)

PRODUCT STORAGE Washing instructions

It can be machine washed at a max temperature of 95°C. Avoid bleaching and ironing. Store the product away from direct sunlight.

MANAGEMENT OF DAMAGED DEVICES

If the device is damaged, notify the manufacturer immediately.

SIDE EFFECTS AND SERIOUS ACCIDENTS

To avoid side effects, the indications for use should be followed carefully. It is a good practice to control your skin daily

to verify the presence of inflammations or flaking. Removal of the compression stocking is recommended in case of skin discolouration, numbness or persistent pain. In the event of a serious accident occurring in connection with the device, the competent authorities of the member state should be notified.

COMPOSITION

- VARISAN® A.T.E.® open toe: polyamide, 76%, élastane 24%. Latex free.
- VARISAN® A.T.E.® with inspection hole: polyamide 77%, élastane 23%. Latex free.
- VARISAN® A.T.E.® with inspection hole COTTON: polyamide 70%, élastane 15%, cotton 15%. Latex free.
- VARISAN® A.T.E.® two-heights, with inspection hole: polyamide 82%, élastane 18%. Latex free.
- VARISAN® RECOVER24: polyamide 62%, élastane 38%. Latex free.
- VARISAN® RECOVER35: polyamide 52%, élastane 48%. Latex free.

VALIDITY

Compression stockings can be worn both during the day and at night, always referring, however, to the instructions of the attending physician and taking care to change them if damaged. They can undergo consecutive washing cycles for a maximum of 10 times. The materials and technologies applied for their technical construction suggest a 6-month period of compression effectiveness for the stockings from their first use. The need for a possible replacement in periods prior to said deadline may arise, if the instructions indicated in this leaflet are not observed.

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CIZETA MEDICALI S.p.A.